

Document information					
ID:	P001-7	Rev:	2	Author:	Tomas von Peltzer
Product nr:				Product name:	ImageBank
Title:	Change Log				
Reviewed by:	Mattias Ekberg	Date:	2014-10-22	Sign:	
Approved by:	Ingvar Falconer	Date:	2014-10-22	Sign:	

## Change Log

This document contains the change log for the ImageBank product, in the form of Release Notes for the official releases of the product.

### Version 2.3:

#### New features

- Support for asterisk i SSN (for journalssystem Alma)
- Date in 'compare' and 'full screen' functions
- Export with date on image name
- Move date into the image instead of outside
- Rectangle Zoom-function
- Archive functionality for images and inform the user when they recover images from archive that they will use the archived images
- Macro-button in full screen view uses auto filter configurations.
- Customize toolbar.
- Possibility to add twain devices to toolbar.
- Removed SSN validation in Patient overview.
- Added plugins to installation.

#### Bug fixes

- Fixed crash when opening all images with F6, then double-click
- Fixed formatting of Danish texts
- Fix of image export from search-view
- Fix of bug when using 'Show patient information' with Danish SSN
- Fix of bug when images imports to wrong place.
- Empty tab in main view
- SSN script in database uses SSN instead of pmid
- Added Danish text in patient card control
- Fix of incorrectly cut image I comparison mode
- Right click on status button in compare view does nothing, removed possibility.
- Fixed Max/Min error message in auto filter when creating new.
- Fixed arithmetic overflow error when moving images with calibration values.
- New possibility to create auto filter for scanner plugis

**Version 2.2.1 (build 28047)**

## New features

- Connect to database wizard from imagebank login with same features as IBSysAdmin. "Create demo license" button added.
- Updated graphics components from atalasoftware.
- Updated GUI components from Telerik.
- Updated installer GUI and additional feature "Demo plugin"

**Version 2.2.1 (CE version):**

## New features

- Detection of an image that has a different scale and different color depth forces re-calibration before measuring.
- If an import fails from an external device a message box with the path to the temporary directory is shown for possible recovery of the image.

## Changes

- Auto filter function needs sysadmin privileges to be accessed.
- Soft-delete implemented on deleting images. Image-file remains in file structure to prevent unintentional loss of data.

**Version 2.2:**

## New features

- Graphics component updates.
- ImageBank now uses .NET 4 framework
- Labeling information in about tab.

## Changes

- New configuration identification "Fix" in version. Now uses Major.Minor.Fix.Build versioning.

## Removed features

- Integrated video capture (this will later be release as a plugin instead)

**Version 2.1**

## 2.1.4701

## New features

- Save calibration option when making a calibration.
- Removed "Is this an X-Ray image" confirm dialog, gray scale images gets automatically labeled as X-Ray.

## 2.1.4527

## New features

- Twain can now use Auto Filter functionality

## 2.1.4355

## New features

- Import function for Auto Filter generated xml to Auto Filter

## 2.1.4220

## New features

- New setting, TWAIN loop function.

## 2.1.3996

New features

- Auto Filters - Shows value on sharpen

Bug Fixes.

- Cut/Paste between studies
- Drag drop (Copy) from "Image Overview" (F6) to study.